**Moulin de laboratoire CD 2**  
**Laboratory Mill CD 2**  
**Molino de laboratorio CD 2**

**Mode d'emploi • User's manual**

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02/2017

<b>Mode d'emploi</b> .....	<b>5</b>
1. Mise en service .....	11
2. Consignes de sécurité .....	12
3. Mouture d'essai .....	13
4. Entretien de l'appareil .....	16
Annexe 1 : Préparation des blés .....	17
Annexe 2 : Tableau d'essai .....	19
<b>Instruction manual</b> .....	<b>21</b>
1. Setup .....	27
2. Safety instructions .....	28
3. Test milling .....	29
4. Maintenance .....	32
Appendix 1 : Conditioning the wheat .....	33
Appendix 2 : Model grinding report .....	35
<b>Modo de empleo</b> .....	<b>37</b>
1. Puesta en servicio .....	43
2. Consignas de seguridad .....	44
3. Moltura de prueba .....	45
4. Mantenimiento del aparato .....	48
Anexo 3 : Preparación de los trigos .....	49
Anexo 4 : Tabla de prueba .....	51

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**GB** **ENGLISH**

### **IMPORTANT**

The technical elements that constitute this manual (text and illustrations) are not contractual, their only target being to bring assistance for using the mill CD 2.

Use, duplication or disclosure of subject data, for any purpose other than relating to the use or servicing of the equipment, is strictly prohibited without the written authorisation of Chopin Technologies.

Laboratory Mill CD 2

## DECLARATION OF CONFORMITY

We, CHOPIN Technologies, address as below,  
declare under our sole responsibility that the following Apparatus :

### ***Laboratory Mill CD2***

When installed and used in accordance with the instructions in the Product Manual,  
is in conformity with the following standards :

- the EEC directive "Electromagnetic compatibility" 2014/30/UE,
- the EEC directive "Low Voltage" 2014/35/UE.

Emmanuel LECOMTE  
Conformance officer  
CHOPIN Technologies

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Laboratory Mill CD 2

## Laboratory Mill CD 2 - Durum wheat - For semolina production

The equipment and its user's guide conform to the EEC directives.

### ■ Warnings

- Use of the equipment for other purposes than those anticipated by the Chopin Technologies is strictly forbidden.
- Thoroughly follow the warnings and instructions given both in the documentation and on the actual equipment.

### ■ Environment

- inside use,
- operating temperature: 10°C to 45°C,
- hygrometry: 10% to 90%,
- power voltage variations not exceeding  $\pm 10\%$  of the nominal voltage,
- pollution grade (following EN 61010-2): 2,
- installation category for transient surges (following EN 61010-2): II.

### ■ Characteristics

- 230/400 Vac tri-phase 50 or 60 Hz (According to the device configuration, consult the data plate.)

and

- 230 Vac mono-phase 50 or 60 Hz (according to the device configuration, consult the data plate.)
- Power consumption:                   for 230/400 V : 1360 W  
                                                          for 230 V : 70 W.
- Net weight:                                   110 kg.
- Dimensions (mm):                        L 1100 x W 450 x H 900

### ■ Electrical connection

- The electrical network must be fitted with a RDC (Residual Current Device) of sensitivity  $\leq 30$  mA of AC class upstream the bench.
- The 2 Ph + Earth plug can only be plugged into a socket fitted with a protection conductor.

### ■ Fuses

Fuse 5 x 20 F 2A 250V

### ■ Transport

Pack the equipment in its original packaging to fully protect it.

Laboratory Mill CD 2

## 1. Setup

### 1.1 Installation

The machine is generally delivered assembled and ready to use.

❑ In certain instances, the two feed-hoppers may be dismantled due to problems of handling bulk during transport. Carefully remove the two distributor protection shield and mount the two feed-hoppers (the one with the magnetic dumper goes on the right or break side; the other on the left side).

❑ If possible, set mill on a sturdy flat-top table (machine weighs approximately 100 kg) measuring 1.20 m by 0.50 m. The height of the table may vary according to the height of the operator. We recommend an average table height of 0.70 m. The table must be placed in the laboratory in such a way so as to allow the operator to freely move around the mill.

### 1.2 Electrical connection

The two 0.5 HP motors and the micro-motor reducers are connected at the back of the machine:

❑ 3 wires and 1 ground for three-phase 230/400 V AC for power supply to the motor. break roll 0.75 CV - reduction roll: 0.5 CV.

❑ 2 wires + 1 ground terminal for the 230 V single phase micro-gear motor supply.

Unless otherwise specified by the customer, the triphase motors are connected to 400 V.

Once connected, check that the 0.5 HP motors are turning in the correct direction: two arrows on each end of the sieves (end plates removed) indicate the correct rotary direction. When the motor on the break roll side is turning in the correct direction, then the reduction roll side motor will also rotate correctly and vice-versa.

❑ In order to connect the CD2 to three-phase 230 V, See "[MOTORS WIRING](#)" page 53.

### 1.3 Preliminary actions prior to start-up

#### ■ Wheat break roll side

- Check that the front and back sight-holes of the grinder are in their proper place.
- Check, by dismantling the sieve end-plate, that the beaters mounted on the shaft turn freely and without coming into hard contact with the sieve mounting (a possibility if there is damage during shipment).
- Check that the receiving pans are carefully placed.

#### ■ Semolina reduction roll side

Same for points referred to above.

## 2. Safety instructions

No contact with the moving parts of the machine is required when using the TEST MILL.

Grinding must be done without manipulating the plugs which close the sight-holes of the machine cylinders and sieve plates.

**Handling any of the moving parts (sight-hole plugs, guide plates) requires stopping machine beforehand.**

**Any eventual repair work involving the central chamber of the machine requires disconnecting the two power supply plugs.**

### 3. Test milling

Milling is carried out using 600 g of durum wheat, which have been cleaned in order to remove all foreign particles. By either dampening or drying, the moisture capacity of the grains is fixed to 16.5 or 17.5 % (see below).

Significant semolina production requires three following operations:

#### 3.1 Conditioning the wheat

##### ■ Cleaning

The grain must be carefully cleaned and all foreign grains (cereals other than wheat and fodder plants), stones and metallic particles removed.

Cleaning is done either by hand (spreading out the sample on a 0.50 x 0.50 m surface) or by an automatic cleaner such the Chopin Technologies sample cleaner "Quatuor" or "NSP".

##### ■ Conditioning to obtain correct moisture level

□ The amount and the distribution of moisture in the grain are extremely important in order to obtain good milling. Thus, it is necessary to normalize moisture content at a fixed level by dampening, and to try to obtain conditions favouring the better yield of semolina.

□ The simplest method is the so-called "cold tempering" technic.

Time of tempering can vary as a function of aims or used technics. For example:

- I.N.R.A. (Institute National of French Agronomical Research)

First 48 hours tempering period to increase moisture to 14.5 %  
second dampening to increase to 17.5 % with 25 min. tempering period.

- Dampening to increase to 16 % with 4 hours tempering period followed by a second dampening of 0.5 % 10 minutes before milling.

- First dampening to increase to 16.5 % with 3-hours tempering period followed by a second dampening to 17.5 % 30 minutes before milling.

The table in appendix 1 shows the amount of water in ml to add to 500 g of grain varying in moisture content in order to bring the water content to 16.5 or 17.5 %.

□ The most common way of preparing wheat in order to obtain the adequate moisture content is to put a sample in 2-litre air-tight flash and shake manually for several minutes in order to completely distribute the water previously added with a graduated tube. This method which has not been standardized in detail appears much too rough when we think of the importance of this preliminary test milling phase.

The addition of water is generally done in a most imprecise way. In particular, the water is unequally distributed over the grain surface. As test milling is an entire process, we felt that it was necessary to complete the. Milling equipment with a fitting designed for:

- improving and standardizing the conditions for preparing the wheat for grinding.
- economy: both in terms of time and actual operation which till now has been performed manually.

With this in mind, a rotary mixer MR 2 L CHOPIN with a water-adding device and a homogenizing screw have been created (refer to the brochure of the machine).

### 3.2 The actual milling

■ Weigh out 500 g of grain in the flask which was used during the preparatory stage. The yields can thus be easily calculated multiplying by 0.2. They correspond to yields obtained in industrial milling with respect to wheat ground at B1.

■ The 500 g of grain are poured through the spout into the feed-hopper at right. The spout which is removable has on the bottom edge a magnet for catching metallic particle which were not removed during the cleaning process.

Before each test, check the spout. If metallic particles are present, remove the spout by unscrewing the star-knob and clean off with a brush.

#### ■ Break

□ Start the grinding by turning the right-hand electric starting switch (clockwise) which simultaneously controls the distributor micro-motor and the grinder and sieve motor.

Make sure that the two electrical connectors are live: starting the distributor without simultaneously starting the grinder will cause jamming.

If this happens,

- Start the machine in the opposite direction.
- Clear away the products with a brush through the two sight-holes.
- Repeat this operation several times until the pulley turns freely.
- Replace the cover and start up again.

□ Grinding time varies with the size of the grains (or the weight of 1,000 grains). This takes about 3 to 3.5 minutes.

After grinding 500 g wait approximately 1.5 min. so that the sieving

is completed. Then, empty the sieve.

□ Two main fractions are obtained:

- at the end, coarse and fine brans,
- coarse semolina in the central collecting-pan.

#### ■ Reduction

□ Pour the semolina into the left feed-hopper of the reduction system (with a flour scoop).

□ Start the motor by turning on the left-hand switch (clockwise) which controls the distributor-reduction roll and sieve unit. Reduction time obviously varies, depending on the amount of semolina extracted during the break stage (around 3/5 minutes).

When the semolina has been reduced, wait approximately 1.5 minutes for the sieving to be completed and for the sieve to empty itself.

□ Three fractions are obtained:

- at the far end, medium coarse semolina,
- in the left collecting-pan, reduction flour,
- in the right collecting pan, fine semolina

■ The time for one grinding of local wheat should not exceed 10 minutes.


The total yield of the milling must be better than 98 % (including moisture loss), the repeatability of the extraction ratio with the same wheat conditioned under identical conditions being  $\pm 1$  %.

Appendix 2 shows a model grinding report.

## 4. Maintenance

The design of the machine avoids complicated maintenance problems.

Simply check out the following points:

■ Pour several drops of vaseline oil on the distributor roller bearings (4 points every 6 months).

■ Regularly check the condition of the sieve fittings. The sieve materials (stainless steel) must not be stressed and must be able to vibrate at the same rotary speed as the sieve shaft. Changing the sieve fittings is done by replacing the complete sieve device.

■ The position of the final beaters has been set during the machine tests. It is strongly recommended to modify this setting only with caution and consideration.

## Appendix 1 : Conditioning the wheat

**Example: Calculation for the quantity of water to add to the 16.5% level**  
(for 17.5% level, follow' same line).

Making a run with wheat having a moisture content A %  
(A % < 16.5 %), for example 10 %, to 16.5 % which is the moisture  
content for the milling test, does not mean that it will be necessary  
to add:

$$16.5 \% - A \%$$

For example,  $16.5 \% - 10 = 6.5 \%$  moisture, or 6.5 ml per 100 g of  
wheat

■ In fact, the composition of the product is written:

$$\text{Dry matter (\%)} + \text{Initial moisture content A (\%)} = 100$$

If one wishes to increase the moisture content % of the product to  
16.5 %, the amount of water to add must take into account the  
amount of dry matter to dampen in order to reach the correct final  
moisture content.

The amount of water is given by the following equation:

$$\frac{A + X}{MS + A + X} = \frac{16,5}{100}$$

or

$$X = \frac{16,5 \times MS}{83,5} - A$$

or again

$$X = (16,5 - A) \times 1,2$$

■ For example: conditioning wheat from 10 to 16.5 % moisture  
content

Quantity of water to add per 100 g

$$X = \frac{16,5 \times 90}{83,5} - 10$$

$$X = \frac{1485}{83,5} - 10 = 17,8 - 10 = 7,8 \text{ ml}$$

or

$$X = (16,5 - 10) \times 1,2 = 6,5 \times 1,2 = 7,8 \text{ ml}$$

**Water adding table in ml  
to condition 500 g of wheat**

Initial moisture content of wheat (%)	Final moisture of wheat	
	16,5%	17,5%
10	39	45
10,5	36	42
11	33	39
11,5	30	36
12	27	33
12,5	24	30
13	21	27
13,5	18	24
14	15	21
14,5	12	18
15	9	15
15,5	6	12
16	3	9
16,5	0	6
17	-	3

For intermediate moisture content, interpolate the above values:

□ Add +0.6 ml per 0.1 %

or use the simplified equation:

□  $X_{ml} = (16.5 - A) \times 6$

(A = initial moisture content %)

## Appendix 2 : Model grinding report

Mill Test No:			
Wheat reference: .....		Moisture before tempering: .....%	
Moisture before milling: .....%			
Quantity milled: ..... g		Tempering time: ..... h	
	g	%	Notes
<b>Break</b>			
Overtails of sieve No 25			
Coarse semolina			
<b>Reduction</b>			
Medium coarse semolina			} Semolina extraction
Fine semolina			
Flour			
<b>Total yield</b>			

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Unless otherwise specified by the customer, the three-phase motors are connected to 400 Vac (Y wiring).

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

❑ In order to connect the CD 2 to three-phase 230 V ( $\Delta$  wiring), proceed as follows:

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- Unplug the CD 2.  
- Remove the higher panel (A) and the front panel (B).

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
- Remove the cover of the connector block of each motor (C).

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
- For each motor, disconnect the power supply cable (D) and remove the 3 connectors (E).

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- Position the 3 connectors as shown on the pictures below and plug the motor

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