EM10

Etuve multicellulaire • Multi compartment oven
• Cámara multicelular

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IMPORTANT

The technical elements that constitute this manual (text and illustrations) are not contractual, their only target being to bring assistance for using the EM10 Multi compartment oven.

Use, duplication or disclosure of subject data, for any purpose other than relating to the use or servicing of the equipment, is strictly prohibited without the written authorisation of Chopin Technologies.

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USER'S MANUAL

DECLARATION OF CONFORMITY CE

We, CHOPIN Technologies, address as below,
declare under our sole responsibility that the following Apparatus :

EM10

When installed and used in accordance with the instructions in the Product Manual,
is in conformity with the following standards :

- the EEC directive "Electromagnetic compatibility" 2014/30/UE,
- the EEC directive "Low Voltage" 2014/35/UE.

Emmanuel LECOMTE
Conformance officer
CHOPIN Technologies

Date :
12/2016

METHODS AND EQUIPMENT FOR THE CONTROL OF CHARACTERISTICS OF CEREALS AND THEIR DERIVATIVES
METHODES ET EQUIPEMENTS POUR LA MAITRISE DES CARACTERISTIQUES DES CEREALES ET DE LEURS DERIVES

Customer Support / SAV Tel + 33 1 41 47 50 33 - service@chopin.fr
CHOPIN TECHNOLOGIES – KPM Analytics Group - 20 avenue Marcellin Berthelot - 92396 Villeneuve-la-Garenne Cedex, France
SAS au capital de 6 478 694 € - RCS Nanterre - Siret 403 156 441 00020 - APE 2651B - TVA FR 03 403 156 441


1. Introduction

■ The EM10

- The EM10 Multi compartment oven is designed to determine the water content of wheat and semolina (and any other cereal).
- The EM 10 is a referenced oven.

■ Principle

- Each compartment of the EM10 is designed to enable strict temperature transmission to any product introduced.
- A **natural** ventilation created by an air flow (intake through door inlets and exhaust through ventilation chimneys) provides for fast elimination of water steam emitted by the drying of samples.

■ Standards

- Meet the following standards:
 - AACC 44-20.01 ; AACC 44-15.02
 - FT WG 0008
 - NF ISO 712 ; NF V03 707 ; NF V03 708
 - IRAM 1580-1
 - ICC 110/1

EM10 Multi compartment oven

■ **General characteristics:**

□ The EM10 consists of 10 compartments with capacity for up to 20 samples.

□ A regulator allows the temperature to vary within $\pm 0.1^{\circ}\text{C}$.

□ Temperature range: from 30°C to 150°C .

□ Power supply: 230 Vca - 50 /60 Hz - 450 w

□ Noise level: < 70 dB

□ Cells dimensions:

Width x Length x Height: 75 x 160 x 40 mm

□ Multi compartment oven dimensions:

Width x Length x Height: 280 x 280 x 580 mm

□ Weight:

- Net / Gross: 14 / 22 kg

■ **Fuse**

Fuse 5x20 T 2A 250V

■ **Environment**

□ Inside use

□ Storage temperature: $- 20$ to $+ 60^{\circ}\text{C}$

□ Operating temperature: 2 to 45°C

□ Hygrometry: $\text{HR} \leq 85\%$ at 40°C

□ Power supply tolerance: $< \pm 10\%$

□ Pollution level according to EN 61010: 2

□ Installation class according to EN 61010: II
(for transitional over voltage)

2. Installation and setup

■ Unpacking



□ The apparatus is delivered in a special packing, that we advise you to keep for future use (transfer of the Multi compartment oven from one site to another one, periodic service...).

This packing contains the EM10, the power cord and the instruction manual.

□ In its delivery packing, the EM10 is maintained in a cardboard equipped with polystyrene foam for preventing shock damage.

To take out the EM10 of its packing, lift up the cardboard by holding its handles.

Put it on a flat surface. Open the vertical sides of the cardboard and free the EM10 from polystyrene foam.


EM10 Multi compartment oven

■ **Composition**

□ The EM10 (1) is delivered with the following equipment:


- 10 basket holders (3);
- baskets (2) (provided separately ; ref EM 55(451798));
- a hook (4);
- a power cord.

■ **Installation**

□ Put the EM10 on a horizontal and stable surface.

□ Avoid placing the EM10 behind a glass window, mainly if this window is exposed to sunshine for several hours a day.

EM10 Multi compartment oven


■ **Electrical connection**

- ❑ The EM10 is an electric instrument, which requires an earth connection. Above all verify the condition of this connection.
- ❑ In addition, the electrical network must incorporate, upline from the unit, a differential circuit breaker with a sensitivity ≤ 30 mA - class AC.
- ❑ The connection to the mains is realized by the means of the power cord delivered with the EM10.
- ❑ The mains plug is on the rear side of the instrument.

■ **Power ON / OFF**

- ❑ The EM10 is not fitted with a ON / OFF switch.
- ❑ The EM10 powering ON / OFF is performed by connecting / disconnecting the apparatus.

3. Use


■ Precautions

During the EM10 operation, stay away from the hot areas.


The EM10 operation requires a minimum of precautions.

To open the Multi compartment oven doors, use the handles.

Do not block the ventilation chimneys.


Do not insert any object in the ventilation holes.

EM10 Multi compartment oven

■ **Setpoint temperature setting**

- The temperature depends on the protocol used.
- The factory preset temperature is 130 °C.
- To **increase** or **decrease** the setpoint, use the regulator keys.


EM10 Multi compartment oven

■ Baskets introduction

□ To introduce baskets in the Multi compartment oven, use the hook and the basket holders.


□ To obtain an homogeneous ventilation through the Multi compartment oven, leave the basket holders (even if empty) in their compartments.

EM10 Multi compartment oven

■ **Baskets storage**

- The EM10 is fitted in its bottom part with a drawer enabling to store up to 20 baskets.
- They have to be arranged as shown below

EM10 Multi compartment oven

4. Servicing and maintenance

This apparatus does not require any particular maintenance.

■ Cleaning

- Once a year, clean the apparatus using a slightly damp cloth.
- Clean the inside of the cells with compressed air.

■ Temperature check

- Each time the temperature setpoint is changed, and at least every six months, check the measured temperature at the level of the superior central screw with a precision thermometer ± 1 °C.

- If there is a temperature difference during this operation, contact your CHOPIN Technologies representative.

■ Fuse replacement

If the apparatus does not switch on and if the power supply has been checked, replace the EM10 fuse.

- Disconnect the apparatus and replace its fuse by using a small screwdriver.
- The fuse is located at the back of the apparatus behind a cover.
- Fuse characteristics: 5x20 T 2A 250V.

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